

### Experiment 1 27 August 2019 Making Measurements and the Density of Pennies

Isn't this

# Objectives: Learn to make the best possible laboratory measurements.

Today we will become familiar with some critical lab skills like picking the best glassware for measurements and making good measurements. We will talk about significant figures, and rounding.

We will collect data and use it to obtain a meaningful result – the density of a metal.

### **Overview:**

- 1. How to use these pre-lab presentations.
- 2. Making good measurements
- 3. Precision, accuracy, significant figures
- 4. Procedure for today (summary of what we do today)
- 5. Density of pennies
  - Your lab report

#### Format of each experiment in the lab manual:

- 1. Background (or Introduction)
- 2. Procedure
- 3. Calculations, Analysis, etc.

Before we get started with Experiment 1, let's buzz a bit about the lab manual and these presentations. Each experiment in the lab manual has three parts as you see in this box. These pre-lab presentations do the same. We try to provide a simple, more concise alternative to the first section (Background). If this works for you, you can skip the Background section in the lab manual.

Buzz. Buzz.

Most of the time, we will closely follow the Procedure from the lab manual. Today... not so much. (See Slide 26.)

And finally, we take you through the calculations and you should follow along with what we are doing.

Each week, we introduce you to the experiment and this should help you write the **Introduction** before coming to lab.

Slide 2 always lists the objectives. These are the talking points for your Introduction.

Uh-huh.

Info for Introduction Looking back to Slide 2 – you can see that we should mention something about making good measurements, precision, accuracy and significant figures and finally density. As you flip through this presentation, look for these bright yellow signs for more info for your Introduction.

4

Introductions usually run about a half-page. They should be in your own words and cover all pertinent topics. Here is an Introduction for Experiment 1. Yours could look somewhat similar but in your own words. Notice how all of the "talking points" were addressed, but only as a preview?

#### Introduction.

This week in lab we will use various glassware, balances and rulers in order to compare how accurately each one makes measurements. We will practice making good measurements and will learn how many significant figures each measurement gives. We will practice reading one more number in between two lines on a ruler or with graduated glassware. We will encounter examples of precision and accuracy with our measurements. In our calculations, we will review how to use significant figures correctly.

We will determine the density of pennies using a s and volume measurements that we will plot so that the slop the density.

These presentations introduce you to the experiment, so you can usually skip the Background/Introduction section from the lab manual and go straight to the Procedure (once we are in lab). These presentations are tailored to go with Dr Mattson's lectures, using the same vocabulary and approaches you've seen in lecture.

The pre-lab presentations also provide an overview of the procedure followed each week, however, they lack the detail that the lab manual gives. So we generally follow the **Procedure** section in the lab manual carefully. Usually, word for word. (But as I mentioned, today is different.)

These pre-lab presentations always walk you through the calculations in a way that is easy to follow. (The lab manual expects you will refer back to the lab manual's introduction for this.)

Bring your laptop and have this presentation available during lab. So when it's time to perform the calculations, you can refer to the presentation. Not every lab will have calculations, but they all have some sort of analysis of how things went.

Let's go!

So this week we are making good measurements?

What about great measurements?

Good (great) measurements start with the proper glassware – or equipment such as balances.

The beaker (that you have your foot in) is not very accurate. The volumetric flask is very accurate – exactly 100.0 mL. Hint: The word volumetric in the name says it's accurate.

See that yellow sign in the corner? That's a hint that there is info on this slide for your introduction. Perhaps something about accuracy and glassware?

8











Real water forms a meniscus – we read from the bottom of the meniscus. ...Aaaah... between 52 and 53... Closer to 53... Sort of 8/10 closer to 53... so... I get 52.8 mL

Photo credit: http://www.austincc.edu/biocr/1406/labm/ex2/images/graduated\_cylinder.jpg



**PRECISION** has two definitions. In our textbook, precision refers to how well a set of independent measurements agree with one another.

In lab, **PRECISION** also refers to how many significant figures the measurement gives. So it has two definitions.

We could measure the mass of a penny multiple times, for example, using various pocket balances that read to the hundredths place.

These data look pretty good. They are all similar so that suggests precision of the first kind.

## The mass of a 2017 penny

Trial 12.76 gTrial 22.74 gTrial 32.73 gTrial 42.74 g

Average 2.743 g

LIBERT

We certainly expect uncertainty in the last decimal place.

16

| The mass of a   |           |  |  |  |  |  |  |
|-----------------|-----------|--|--|--|--|--|--|
| 2017 per        | iny using |  |  |  |  |  |  |
| pocket balances |           |  |  |  |  |  |  |
| Trial 1         | 2.76 g    |  |  |  |  |  |  |
| Trial 2         | 2.74 g    |  |  |  |  |  |  |
| Trial 3         | 2.73 g    |  |  |  |  |  |  |
| Trial 4         | 2.74 g    |  |  |  |  |  |  |
|                 |           |  |  |  |  |  |  |

Average 2.743 g



Using analytical balances we get one or two extra places past the decimal point – better precision (of the second kind)

> The mass of the same 2017 penny using different analytical balances

> > Trial 12.763 gTrial 22.751 gTrial 32.729 gTrial 42.744 g

Average 2.7468 g

**ACCURACY** refers to how well the values compare to the actual true value, assuming it is known.

Ooops

If the four trials were made by four people using four different balances, and the data were consistent, we would be pretty certain of the accuracy.

> Suppose four people used the same balance – and maybe it is not working right. Or suppose one person measured a volume four times but was doing it wrong. Accuracy is not certain.

 That brings us to significant figures
and the rules for adding, subtracting, multiplying, and dividing them.

How many significant figures are present in each of these values?

Ooo ooo! Ask me! I'm a numbers bird. 18.99 0.0042 0.01610 12.002 47000 100.

How did you do?

18.99 (FOUR) 0.0042 (TWO) 0.01610 (FOUR) 12.002 (FIVE) 47000 (TWO) 100. (THREE)

Info or Introduction

Here are a two examples of conversions that are not exact - both have five sig figs:

1 L = 1.0567 qt 1 pound = 453.59 g Moving along, numbers that are **defined** have an infinite number of significant figures because they are exact. For example, 1 dozen = 12, or 1 m = 100 cm. And my personal favorite: 1 mL = 1 cm<sup>3</sup>.

> Here are two bizarre exact conversions (by definitions): 1 inch = 2.54 cm exactly 1 cal = 4.184 J exactly

The rule for adding or subtracting numbers is trickier than the rule for multiplying or dividing... We'll start with tricky.

Here's how: Line up the decimal points of the numbers to be added or subtracted. This involves writing the numbers down all lined up by decimal point. Then add them up using your calculator (or your head). Chop off the answer to match the number that has the fewest places past the decimal. EXAMPLE 1 8.65 4 <u>+ 94.03</u> 102.68

See? You can easily get more (or fewer) significant figures than you started with!

Example 2

2.1

+ 7.422

9.5

9.522

| Evample 3            | Evample 1  |
|----------------------|------------|
|                      | Litample 4 |
| 0.00 <mark>77</mark> | 0.0027     |
| + 1.04               | + 1.040    |
| 1.04 <mark>77</mark> | 1.0427     |
| 1.05                 | 1.043      |

Here are some more examples. See how we rounded in Examples 3 and 4? Also, listen up! This is important! Keep your full answer in your calculator until the very end – just keep track of how many figures your answer should have. This prevents propagating rounding errors. In Example 2, keep 9.522 for the next calculation, but remember you get only two sig figs from it in the final answer..

*The rule for multiplying and dividing... Well, we all know that rule...* 

If one number has five significant figures and the other number has only three, the product or quotient has only three – the same as the least of the participant numbers.

Example 1. A density problem: Suppose 2.6943 g has a volume of 1.37 mL. What is its density?

density = mass/volume = 2.6943 g / 1.37 mL = 1.966642336 g/mL → 1.97 g/mL

Again, keep the whole number in your calculator until the end of the calculations Just remember there are supposed to be three significant figures.

# Example 2. What is the average mass of the penny from the pocket balance measurements?

Here we use both rules: We add up the four masses using the adding rule and then divide by 4 (an exact number this time) using the dividing rule. The mass of a 2017 penny using pocket balance Trial 1 2 76 g

|         | 2.70 g        |
|---------|---------------|
| Trial 2 | <b>2.74</b> g |
| Trial 3 | <b>2.73</b> g |
| Trial 4 | 2.74 g        |

Trial 12.76 gTrial 22.74 gTrial 32.73 gTrial 42.74 g

Sum = 10.97 g

**4 significant figures** 

DIVIDE BY 4 (EXACT) Average 2.743 g

LIBERT

Info for calculations



### 4. Procedure: What we will do today

Now it's time for you to do the experiment. The TAs will help you with the lab safety part of the Procedure. You do not need to write anything in your lab notebook.

The TAs will also help you with the next two parts, Glassware and Equipment and practicing taking measurements. Instead of the lab manual procedure for the density of pennies, follow the procedure on the back of the cover sheet and overviewed on the next four slides.

> Do the last part of the lab manual procedure on precision & accuracy if you have enough time. It will not be graded.

### 5. Procedure for the density experiment. (Details given on the back of the cover sheet.)



Your experiment will look like this sequence as it goes along. Start with five piles of 10 pennies each and weigh each pile to the nearest 0.01 g. Fill a graduated cylinder with about 40 mL water. Read the volume to the nearest 0.1 mL. Add the first pile of pennies and read the volume again. Record the mass of pennies just added and the new volume. Repeat until all five piles have been added.

### 5. Procedure for today (our density experiment)

 Follow the instructions on the back of your cover sheet to create an Excel spreadsheet with your data.

| 1 |       | В               | С              | D            | E               | F                   | G                |  |
|---|-------|-----------------|----------------|--------------|-----------------|---------------------|------------------|--|
| 1 | Trial | Number of Coins | Initial Volume | Final Volume | Volume of Coins | Mass of coins added | Total mass coins |  |
| 2 | 1     | 10              | 40.1           | 44.0         | 3.9             | 27.0                | 27.02            |  |
| 3 | 2     | 20              | 40.1           | 47.7         | 7.6             | 27.6                | 54.60            |  |
|   | 3     | 30              | 40.1           | 51.0         | 10.9            | 27.0                | 81.62            |  |
| 5 | 4     | 40              | 40.1           | 54.9         | 14.8            | 27.1                | 108.74           |  |
| 6 | 5     | 50              | 40.1           | 58.1         | 18.0            | 27.1                | 135.83           |  |
| 7 |       |                 |                |              |                 |                     |                  |  |

Learn to write formulas, learn to copy and paste formulas and how to format your numbers. Help each other and ask Dr. Mattson or the TAs if you have questions.

You and your lab partner can work together on this, but print a copy for each of you.

<u>)</u>8

### 5. Procedure for today (our density experiment)





### 6. What your lab report should include.

- First, the cover sheet with TA initials.
- You now know how to write an Introduction (Slide 5). We will discuss conclusions and sources of error more in the pre-lab to Experiment 3.
- (3)Next, attach the copy pages from your lab notebook (trim off rough edges)
- (4)Last, the printed Excel x-y scatter graph and data from our density of pennies experiment. Each person needs a graph.
- (5) All of the above stapled together.
- (6)Turned in lab report (each individual) to correct folder by start of class tomorrow. Or... turn it in right after lab
  - today. Late labs may not be graded see the syllabus Next week: We will need to

dress for a mess.

Stick people inspired by xkcd cartoons by Randall Munroe (www.xkcd.com)

2)

Chem Lab with the Stick People and Bird was created and produced by Dr. Bruce Mattson, Creighton Chemistry. Enjoy it and share it if you wish.

31

Woo-hoo!

We did it!

First lab!